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β-Laktamase für Staphylokokken EUCAST/CLSI

Staphylococcus spp.


EUCAST Clinical Breakpoint Table v. 4.0, valid from 2014-01-01

<p>Disk diffusion (EUCAST standardised disk diffusion method) Medium: Mueller-Hinton agar Inoculum: McFarland 0.5 Incubation: Air, 35±1°C, 18±2h Reading: Read zone edges as the point showing no growth viewed from the back of the plate against a dark background illuminated with reflected light. Quality control: <i>Staphylococcus aureus</i> ATCC 29213</p>
--

Penicillins ¹	MIC breakpoint (mg/L)		Disk content (µg)	Zone diameter breakpoint (mm)		Notes Numbers for comments on MIC breakpoints Letters for comments on disk diffusion
	S ≤	R >		S ≥	R <	
						1/A. Most staphylococci are penicillinase producers, which are resistant to benzylpenicillin, phenoxymethylpenicillin, ampicillin, amoxicillin, piperacillin and ticarcillin. Isolates negative for penicillinase and susceptible to methicillin can be reported susceptible to these agents. Isolates positive for penicillinase and methicillin susceptible are susceptible to beta-lactamase inhibitor combinations and isoxazolypenicillins (oxacillin, cloxacillin, dicloxacillin and flucloxacillin). Methicillin resistant isolates are, with few exceptions, resistant to all beta-lactam agents.
Benzylpenicillin, <i>S. aureus</i>	0.12 ¹	0.12 ^{1,2}	1 unit	26 ^{A,B}	26 ^{A,B}	B. For <i>S. aureus</i> , disk diffusion is more reliable than MIC determination for detection of penicillinase producers, provided the zone diameter is measured AND the zone edge closely inspected (see pictures below). If the zone diameter is <26 mm, then report resistant. If the zone diameter is ≥26 mm AND the zone edge is sharp, then report resistant. If not sharp, then report susceptible and if uncertain, then report resistant. Chromogenic cephalosporin-based beta-lactamase tests do not reliably detect staphylococcal penicillinase.
Benzylpenicillin, <i>S. lugdunensis</i>	0.12 ¹	0.12 ^{1,2}	1 unit	26 ^A	26 ^A	
Benzylpenicillin, Coagulase negative staphylococci	- ³	- ³		Note ^C	Note ^C	3/C. No currently available method can reliably detect penicillinase production in coagulase-negative staphylococci.
Ampicillin, <i>S. saprophyticus</i>	Note ¹	Note ¹	2	18 ^{A,D}	18 ^{A,D}	4/D. Ampicillin susceptible <i>S. saprophyticus</i> are <i>mecA</i> -negative and susceptible to ampicillin, amoxicillin and piperacillin (without or with a beta-lactamase inhibitor).

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


Examples of inhibition zones for *Staphylococcus aureus* with benzylpenicillin.

- a) Fuzzy zone edge and zone diameter ≥ 26 mm. Report susceptible.
b) Sharp zone edge and zone diameter ≥ 26 mm. Report resistant.

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JOURNAL OF CLINICAL MICROBIOLOGY, June 2002, p. 2153–2162
0095-1137/02/\$04.00+0 DOI: 10.1128/JCM.40.6.2153–2162.2002
Copyright © 2002, American Society for Microbiology. All Rights Reserved.

Vol. 40, No. 6

Detection of Plasmid-Mediated AmpC β -Lactamase Genes in Clinical Isolates by Using Multiplex PCR

F. Javier Pérez-Pérez^{1,2} and Nancy D. Hanson^{1*}

Center for Research in Anti-Infectives and Biotechnology, Department of Medical Microbiology and Immunology, School of Medicine, Creighton University, Omaha, Nebraska,¹ and Department of Immunology, Microbiology, and Parasitology, Basque Country University, Vitoria-Gasteiz, Spain²

- ACC *Hafnia alvei*
- FOX *Aeromonas* spp.
- MOX *Aeromonas* spp.
- DHA *Morganella morganii*
- CIT *Citrobacter freundii*
- EBC *Enterobacter cloacae*


FIG. 1. AmpC dendrogram. Sequences were downloaded from the GenBank database, and structural genes were compared, as described in Material and Methods, by using the DNAsis program. Values in blue correspond to the percent similarity between the most distinct member of each cluster and the other members within that cluster. Primer pairs are correlated by the family of genes that they amplify.

2015

1. *Escherichia coli* **SAC Amoxicillin/Clavulansäure 16-18 mm intern.**
2. *Staphylococcus aureus*
3. *Rothia mucilaginosa*
4. *Corynebacterium diphtheriae*
Toxin negativ
5. *Streptococcus mitis*
hoher Penicillin - Resistenz
6. *Staphylococcus saprophyticus*
EUCAST Ampicillin testen
7. *S. aureus* **MRSA**
Beta-Laktamase -
8. *Rhizobium radiobacter*
9. *Salmonella* Enteritidis
10. *Neisseria zoodegmatis* (EF-4b)
11. *Proteus mirabilis* **ESBL - CTX-M,**
gleichzeitig **AmpC - CIT-M**
12. *Hafnia alvei*
13. *Neisseria gonorrhoeae*
DD: *Kingella denitrificans*
14. *Streptococcus pseudoporcinus* DD:
Streptococcus agalactiae
15. *Inquilinus limosus* bei CF
16. *Pantoea agglomerans*
17. *Klebsiella pneumoniae* **ESBL ↑↑**
Mit Porindefekt? → Ertapenem res.
18. *Corynebacterium jeikeium*
19. *Arcanobacterium haemolyticum* –
abgeschwächter CAMP
20. *Kerstersia gyiorum* – NF Beinwunde

Kingella denitrificans in Urethralabstrich

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
Jost G, Bloemberg GV, Hombach M. Improved sensitivity for methicillin resistance detection in coagulase-negative staphylococci by moxalactam antibiotic discs or a cefoxitin investigation zone.

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Isolation of *Streptococcus dysgalactiae* subsp. *equisimilis* possessing Lancefield's group A antigen from two patients

V. Deggim, R. Zbinden

Institute of Medical Microbiology, University of Zurich, 8006 Zurich, Switzerland

Results (1)

Isolate 1 and Isolate 2 formed large colonies and were beta-hemolytic (Figure 1).


Figure 1: Growth of Isolate 1 and Isolate 2 in comparison to a *S. pyogenes* strain on a sheep blood agar plate after overnight incubation at 37°C.


Figure 2: Lancefield's group antigen identification of Isolate 1 and Isolate 2.

	Isolate 1	Isolate 2
Hemolysis	Beta- hemolytic	Beta- hemolytic
Colonies	Large colonies	Large colonies
Lancefield's group antigen	A	A
PYR	negative	negative
Bacitracin susceptibility	resistant	resistant
VP reaction	negative	negative
Vitek 2 GP card identification	<i>S. dysgalactiae</i> subsp <i>equisimilis</i> (excellent identification)	<i>S. dysgalactiae</i> subsp <i>equisimilis</i> (very good identification)

Table 1. Phenotypic and biochemical characteristics of Isolate 1 and Isolate 2.

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
1. *Escherichia coli* **ESBL**
2. *Staphylococcus lugdunensis*
MLS und Methicillin-Resistenz
3. *Burkholderia stabilis*
4. *Clostridium tertium*
5. *Ochrobactrum intermedium*
6. *Staphylococcus saprophyticus*
Ampicillin testen
Cefoxitin Grenzwert 22 mm
7. *Morganella morganii*
Imipenem ist nicht mehr
möglich als empfindlich
8. *Clostridium perfringens*
9. *Candida albicans*
10. *Nocardia farcinica*
11. *Klebsiella oxytoca* **DHA induz. AmpC**
12. *Staphylococcus aureus* **MRSA mecC**
13. *S. aureus* atypisch gelb
14. *S. aureus* atypisch trocken
15. *Candida auris*
16. *Aerococcus urinae*
17. *Kluyvera* sp. **CTX-M ESBL**
18. *Sphingomonas paucimobilis*
19. *Corynebacterium urealyticum*
20. *Streptococcus gallolyticus (bovis)*

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
Staphylococcus spp.

Expert Rules and Expected Phenotypes

For abbreviations and explanations of breakpoints, see the Notes sheet

EUCAST Clinical Breakpoint Tables v. 13.1, valid from 2023-06-29

Aminoglycosides ¹	MIC breakpoints (mg/L)			Disk content (µg)	Zone diameter breakpoints (mm)			Notes
	S ≤	R >	ATU		S ≥	R <	ATU	
Amikacin ² , <i>S. aureus</i>	(16) ¹	(16) ¹		30	(15) ^A	(15) ^A		1/A. For information on how to use breakpoints in brackets, see https://www.eucast.org/eucastguidancedocuments/ . 2. Resistance to amikacin is most reliably determined by testing with kanamycin (MIC >8 mg/L). The corresponding zone diameter for the kanamycin 30 µg disk is R<18 mm for <i>S. aureus</i> and R<22 mm for coagulase-negative staphylococci.
Amikacin ² , Coagulase-negative staphylococci	(16) ¹	(16) ¹		30	(15) ^A	(15) ^A		
Gentamicin, <i>S. aureus</i>	(2) ¹	(2) ¹		10	(18) ^A	(18) ^A		
Gentamicin, Coagulase-negative staphylococci	(2) ¹	(2) ¹		10	(22) ^A	(22) ^A		
Netilmicin	IE	IE			IE	IE		
Tobramycin, <i>S. aureus</i>	(2) ¹	(2) ¹		10	(18) ^A	(18) ^A		
Tobramycin, Coagulase-negative staphylococci	(2) ¹	(2) ¹		10	(20) ^A	(20) ^A		

Glycopeptides and lipoglycopeptides ¹	MIC breakpoints (mg/L)			Disk content (µg)	Zone diameter breakpoints (mm)			Notes
	S ≤	R >	ATU		S ≥	R <	ATU	
Dalbavancin ²	0.125 ^{3,4}	0.125 ³			Note ^A	Note ^A		1. Glycopeptide MICs are method dependent and should be determined by broth microdilution (ISO standard 20776-1). <i>S. aureus</i> with vancomycin MIC values of 2 mg/L are on the border of the wild-type distribution and there may be an impaired clinical response. 2. Resistant isolates are rare or not yet reported. The identification and antimicrobial susceptibility test result on any such isolate must be confirmed and the isolate sent to a reference laboratory. 3. MICs must be determined in the presence of polysorbate-80 (0.002% in the medium for broth dilution methods; agar dilution methods have not been validated). Follow the manufacturers' instructions for commercial systems. 4. <i>S. aureus</i> isolates susceptible to vancomycin can be reported susceptible to dalbavancin and oritavancin. 5. MRSA isolates susceptible to vancomycin can be reported susceptible to telavancin.
Oritavancin ² , <i>S. aureus</i>	0.125 ^{3,4}	0.125 ³			Note ^A	Note ^A		
Teicoplanin ² , <i>S. aureus</i>	2	2			Note ^A	Note ^A		
Teicoplanin, Coagulase-negative staphylococci	4	4			Note ^A	Note ^A		
Telavancin ² , MRSA	0.125 ^{3,5}	0.125 ³			Note ^A	Note ^A		
Vancomycin ² , <i>S. aureus</i>	2	2			Note ^A	Note ^A		
Vancomycin ² , Coagulase-negative staphylococci	4	4			Note ^A	Note ^A		

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Staphylococcus spp.

EUCAST Clinical Breakpoint Tables v. 13.1, valid from 2023-06-29

Expert Rules and Expected Phenotypes

For abbreviations and explanations of breakpoints, see the Notes sheet

Penicillins ¹	MIC breakpoints (mg/L)			Disk content	Zone diameter breakpoints (mm)			Notes
	S ≤	R >	ATU		S ≥	R <	ATU	
Benzylpenicillin, <i>S. aureus</i>	0.125 ¹	0.125 ¹		1 unit	26 ^{A,B}	26 ^{A,B}		Numbered notes relate to general comments and/or MIC breakpoints. Lettered notes relate to the disk diffusion method. 1/A. Most <i>S. aureus</i> are penicillinase producers and some are methicillin resistant. Either mechanism renders them resistant to benzylpenicillin, phenoxymethylpenicillin, ampicillin, amoxicillin, piperacillin and ticarcillin. Isolates that test susceptible to benzylpenicillin and ceftiofur can be reported susceptible to all penicillins. Isolates that test resistant to benzylpenicillin but susceptible to ceftiofur are susceptible to β-lactam β-lactamase inhibitor combinations, the isoxazolympenicillins (oxacillin, cloxacillin, dicloxacillin and flucloxacillin) and nafcillin. For agents given orally, care to achieve sufficient exposure at the site of the infection should be exercised. Isolates that test resistant to ceftiofur are resistant to all penicillins. 2/C. Most staphylococci are penicillinase producers and some are methicillin resistant. Either mechanism renders them resistant to benzylpenicillin, phenoxymethylpenicillin, ampicillin, amoxicillin, piperacillin and ticarcillin. No currently available method can reliably detect penicillinase production in all species of staphylococci but methicillin resistance can be detected with ceftiofur as described. 3/D. Ampicillin susceptible <i>S. saprophyticus</i> are <i>mecA</i> -negative and susceptible to ampicillin, amoxicillin and piperacillin (without or with a beta-lactamase inhibitor). 4. <i>S. aureus</i> , <i>S. lugdunensis</i> and <i>S. saprophyticus</i> with oxacillin MIC values >2 mg/L are mostly methicillin resistant due to the presence of the <i>mecA</i> or <i>mecC</i> gene. Occasionally oxacillin MIC values are high in <i>S. aureus</i> in absence of <i>mec</i> -gene mediated resistance. These isolates have been called BORSA (borderline oxacillin resistant <i>S. aureus</i>). EUCAST does not recommend systematic screening for BORSA. For coagulase-negative staphylococci other than <i>S. saprophyticus</i> and <i>S. lugdunensis</i> , the oxacillin MIC in methicillin resistant isolates is >0.25 mg/L. B. For <i>S. aureus</i> , disk diffusion is more reliable than MIC determination for detection of penicillinase producers, provided the zone diameter is measured AND the zone edge for isolates with zone diameters ≥26 mm is closely inspected (see note 1). If the zone diameter is <26 mm, the zone edge is not inspected (see note 2). If the zone diameter is ≥26 mm, the zone edge is inspected like a "beach", then report lactamase tests do not.
Benzylpenicillin, <i>S. lugdunensis</i>	0.125	0.125		1 unit	26	26		
Benzylpenicillin, other staphylococci	Note ²	Note ²			Note ^C	Note ^C		
Ampicillin, <i>S. saprophyticus</i>	Note ^{2,3}	Note ^{2,3}		2	18 ^{C,D}	18 ^{C,D}		
Ampicillin-sulbactam	Note ^{1,2,3}	Note ^{1,2,3}			Note ^{A,C,D}	Note ^{A,C,D}		
Amoxicillin	Note ^{1,2,3}	Note ^{1,2,3}			Note ^{A,C,D}	Note ^{A,C,D}		
Amoxicillin-clavulanic acid	Note ^{1,2,3}	Note ^{1,2,3}			Note ^{A,C,D}	Note ^{A,C,D}		
Piperacillin	Note ^{1,2,3}	Note ^{1,2,3}			Note ^{A,C,D}	Note ^{A,C,D}		
Piperacillin-tazobactam	Note ^{1,2,3}	Note ^{1,2,3}			Note ^{A,C,D}	Note ^{A,C,D}		
Ticarcillin	Note ^{1,2}	Note ^{1,2}			Note ^{A,C}	Note ^{A,C}		
Ticarcillin-clavulanic acid	Note ^{1,2}	Note ^{1,2}			Note ^{A,C}	Note ^{A,C}		
Temocillin	-	-			-	-		
E. For screening for methicillin resistance in <i>S. pseudintermedius</i>, <i>S. schleiferi</i> and <i>S. coagulans</i>.								
Coagulase-negative staphylococci				1	20 ^E	20 ^E		
Oxacillin (screen only), <i>S. pseudintermedius</i> , <i>S. schleiferi</i> and <i>S. coagulans</i>	NA	NA			Note ^A	Note ^A		
Oxacillin ⁴ , other staphylococci	Note ^{1,4}	Note ^{1,4}			Note ^A	Note ^A		

4. *S. aureus*, *S. lugdunensis* and *S. saprophyticus* with oxacillin MIC values >2 mg/L are mostly methicillin resistant due to the presence of the *mecA* or *mecC* gene. Occasionally oxacillin MIC values are high in *S. aureus* in absence of *mec*-gene mediated resistance. These isolates have been called BORSA (borderline oxacillin resistant *S. aureus*). EUCAST does not recommend systematic screening for BORSA. For coagulase-negative staphylococci other than *S. saprophyticus* and *S. lugdunensis*, the oxacillin MIC in methicillin resistant isolates is >0.25 mg/L.

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Staphylococcus spp.

Expert Rules and Expected Phenotypes

For abbreviations and explanations of breakpoints, see the Notes sheet

EUCAST Clinical Breakpoint Tables v. 13.1, valid from 2023-06-29

Cephalosporins ¹	MIC breakpoints (mg/L)			Disk content (µg)	Zone diameter breakpoints (mm)			Notes
	S ≤	R >	ATU		S ≥	R <	ATU	
Cefaclor ²	Note ¹	Note ¹			Note ^A	Note ^A		Numbered notes relate to general comments and/or MIC breakpoints. Lettered notes relate to the disk diffusion method. 1/A. Susceptibility of staphylococci to cephalosporins is inferred from the cefoxitin susceptibility except for cefixime, ceftazidime, ceftazidime-avibactam, ceftibuten and ceftolozane-tazobactam, which do not have breakpoints and should not be used for staphylococcal infections. For agents given orally, care to achieve sufficient exposure at the site of the infection should be exercised. If cefotaxime and ceftriaxone are reported for methicillin-susceptible staphylococci, these should be reported "Susceptible, increased exposure" (I). Some methicillin-resistant <i>S. aureus</i> are susceptible to ceftaroline and ceftobiprole, see Notes 5/D and 7/F. 2. See table of dosages. 3. <i>S. aureus</i> and <i>S. lugdunensis</i> with cefoxitin MIC values >4 mg/L and <i>S. saprophyticus</i> with cefoxitin MIC values >8 mg/L are methicillin resistant, mostly due to the presence of the <i>mecA</i> or <i>mecC</i> gene. Disk diffusion reliably predicts methicillin resistance. 4. For staphylococci other than <i>S. aureus</i> , <i>S. lugdunensis</i> and <i>S. saprophyticus</i> , the cefoxitin MIC is a poorer predictor of methicillin resistance than the disk diffusion test. 5/C. In <i>S. pseudintermedius</i> , <i>S. schleiferi</i> and <i>S. coagulans</i> the cefoxitin disk is less predictive for the detection of methicillin resistance than in other staphylococci. Use the oxacillin 1 µg disk with zone diameter breakpoints S≥20, R<20 mm. 6/D. Methicillin-susceptible isolates can be reported susceptible to ceftaroline without further testing. 7/E. Resistant isolates are rare. 8/F. Methicillin-susceptible isolates can be reported susceptible to ceftobiprole without further testing. B. If coagulase-negative staphylococci are not identified to species level, use zone diameter breakpoints S≥25, R<25 mm, with an ATU of 22-24 mm. For isolates with results inside the ATU: identify species, perform PCR for <i>mecA/mecC</i> or report resistant.
Cefadroxil	Note ¹	Note ¹			Note ^A	Note ^A		
Cefalexin	Note ¹	Note ¹			Note ^A	Note ^A		
Cefazolin	Note ¹	Note ¹			Note ^A	Note ^A		
Cefepime	Note ¹	Note ¹			Note ^A	Note ^A		
Cefiderocol	-	-			-	-		
Cefixime	-	-			-	-		
Cefotaxime ²	Note ¹	Note ¹			Note ^A	Note ^A		
Cefoxitin (screen only), <i>S. aureus</i> and coagulase-negative staphylococci except <i>S. epidermidis</i> and <i>S. lugdunensis</i>	Note ^{3,4}	Note ^{3,4}		30	22 ^{A,B}	22 ^{A,B}		
Cefoxitin (screen only), <i>S. epidermidis</i> and <i>S. lugdunensis</i>	Note ⁴	Note ⁴		30	27 ^{A,B}	27 ^{A,B}	27	
Cefoxitin (screen only), <i>S. pseudintermedius</i> , <i>S. schleiferi</i> and <i>S. coagulans</i>	Note ⁵	Note ⁵			Note ^C	Note ^C		
Cefpodoxime	Note ¹	Note ¹			Note ^A	Note ^A		
Ceftaroline (indications other than pneumonia), <i>S. aureus</i>	1 ⁶	2 ^{6,7}	1	5	20 ^D	17 ^{D,E}	19-20	
Ceftaroline (pneumonia), <i>S. aureus</i>	1 ⁶	1 ⁶	1	5	20 ^D	20 ^D	19-20	
Ceftazidime	-	-			-	-		
Ceftazidime-avibactam	-	-			-	-		
Ceftibuten	-	-			-	-		
Ceftobiprole, <i>S. aureus</i>	2 ⁸	2 ⁸	2	5	17 ^F	17 ^F	16-17	
Ceftolozane-tazobactam	-	-			-	-		
Ceftriaxone ²	Note ¹	Note ¹			Note ^A	Note ^A		
Cefuroxime iv	Note ¹	Note ¹			Note ^A	Note ^A		
Cefuroxime oral	Note ¹	Note ¹			Note ^A	Note ^A		

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Informations needed for the expert laboratories to communicate to Anresis

Isolate Lab ID (identifier) _____
Patient initials (first name, surname) __ __ date of birth (dd/mm/yyyy): _____ gender f m
Name / address of the first laboratory _____

Sampling date If not available, replace with arrival date in the laboratory _____
Type of specimen blood sterile puncture urine stool
 tracheobronchial secretions other unknown
description: _____
Provider Hospital: _____
 ward ICU emergency room outpatient clinic unknown
 Practice (Name, ZIP-code): _____
 Other: _____
Clinical setting: clinical importance: colonization infection unknown
Clinical description: _____ unknown
epidemiological context: hospital acquired community onset unknown
(after 48 hours)
14-day-mortality: yes no unknown
previous travel or stay abroad within last 6 months: country _____
 yes - hospitalized yes - not hosp. no unknown

Results of resistance mechanism determined by the expert laboratory

Species _____
Carbapenemase _____
Conventional results _____
Molecular results _____
Date: _____ Signature: _____

Transfer of the results to Anresis by SAC: Date: _____ Signature: _____

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Penicillins	MIC breakpoints (mg/L)			Disk content (µg)	Zone diameter breakpoints (mm)			Notes
	S ≤	R >	ATU		S ≥	R <	ATU	
Benzympenicillin	-	-		-	-			Numbered notes relate to general comments and/or MIC breakpoints. Lettered notes relate to the disk diffusion method. 1. For information on how to implement the new aminopenicillin breakpoints, see https://www.eucast.org/eucastguidancedocuments/ . 2. For susceptibility testing purposes, the concentration of sulbactam is fixed at 4 mg/L. 3/D. For information on how to use breakpoints in brackets, see https://www.eucast.org/eucastguidancedocuments/ . 4. For susceptibility testing purposes, the concentration of clavulanic acid is fixed at 2 mg/L. 5. For susceptibility testing purposes, the concentration of tazobactam is fixed at 4 mg/L. 6. Agar dilution is the reference method for mecillinam MIC determination. A. Ignore growth that may appear as a thin inner zone on some batches of Mueller-Hinton agars. B. Susceptibility inferred from ampicillin (iv or oral). C. Isolates susceptible to ampicillin (iv or oral) can be reported "susceptible, increased exposure" (I) to "ampicillin oral (infections originating from the urinary tract)". Isolates resistant to ampicillin (iv or oral) can be reported resistant to "ampicillin oral (infections originating from the urinary tract)". E. Infer from ampicillin oral, but the report should explain the meaning of breakpoints in brackets. F. Ignore isolated colonies within the inhibition zone.
Ampicillin iv ¹	8	8		10	14 ^A	14 ^A		
Ampicillin oral (uncomplicated UTI only) ¹	8	8		10	14 ^A	14 ^A		
Ampicillin-sulbactam iv ¹	8 ²	8 ²		10-10	14 ^A	14 ^A		
Ampicillin-sulbactam oral (uncomplicated UTI only) ¹	8 ²	8 ²		10-10	14 ^A	14 ^A		
Amoxicillin iv ¹	8	8		-	Note ^B	Note ^B		
Amoxicillin oral (infections originating from the urinary tract) ¹	0.001	8		-	Note ^C	Note ^C		
Amoxicillin oral (uncomplicated UTI only) ¹	8	8		-	Note ^B	Note ^B		
Amoxicillin oral (other indications) ¹	(8) ³	(8) ³		-	Note ^{D,E}	Note ^{D,E}		
Amoxicillin-clavulanic acid iv ¹	8 ⁴	8 ⁴		20-10	19 ^A	19 ^A	19-20	
Amoxicillin-clavulanic acid oral (infections originating from the urinary tract) ¹	0.001 ⁴	8 ⁴		20-10	50 ^A	19 ^A	19-20	
Amoxicillin-clavulanic acid oral (uncomplicated UTI only) ¹	32 ⁴	32 ⁴		20-10	16 ^A	16 ^A		
Amoxicillin-clavulanic acid oral (other indications) ¹	(8) ^{3,4}	(8) ^{3,4}		20-10	(19) ^{A,D}	(19) ^{A,D}	19-20	

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