



## Survey report MQ 2016-4 U3 Drugs of abuse

### Amphetamine

	Target value	Participant	Conform	Conform%
Architect	Negative	1	1	100
Roche Cobas	Negative	3	3	100
Triage 8	Negative	8	8	100
Syva Rapid Test	Negative	2	2	100
Intex	Negative	5	5	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	12	12	100
Other quantitativ	Negative	13	13	100
Other Rapid tests	Negative	49	47	96
<b>TOTAL</b>		<b>94</b>	<b>92</b>	<b>98</b>

### Methamphetamine

	Target value	Participant	Conform	Conform%
Intex	Negative	2	2	100
Syva Rapid Test	Negative	4	4	100
Triage TOX	Negative	12	12	100
Other Rapid tests	Negative	23	22	96
<b>TOTAL</b>		<b>41</b>	<b>40</b>	<b>98</b>

### Ecstasy pos/neg

	Target value	Participant	Conform	Conform%
Other Rapid tests	Negative	<b>13</b>	<b>13</b>	<b>100</b>

### Barbiturate

	Target value	Participant	Conform	Conform%
Architect	Negative	1	1	100
Cobas	Negative	3	3	100
Triage 8	Negative	8	8	100
Intex	Negative	1	1	100
Triage TOX	Negative	12	12	100
Other quantitativ	Negative	9	9	100
Other Rapid tests	Negative	23	23	100
<b>TOTAL</b>		<b>57</b>	<b>57</b>	<b>100</b>

**Benzodiazepine**

	Target value	Participant	Conform	Conform%
Architect	Negative	1	1	100
Cobas	Negative	4	4	100
Triage 8	Negative	8	8	100
Syva Rapid Test	Negative	1	1	100
Intex	Negative	10	10	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	12	12	100
Other quantitativ	Negative	14	14	100
Other Rapid tests	Negative	56	55	98
<b>TOTAL</b>		<b>107</b>	<b>106</b>	<b>99</b>

**Cannabinoide**

	Target value	Participant	Conform	Conform%
Architect	Negative	1	1	100
Cobas	Negative	6	6	100
Triage 8	Negative	8	8	100
Syva Rapid Test	Negative	7	7	100
Intex	Negative	11	11	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	12	12	100
Other quantitativ	Negative	15	15	100
Other Rapid tests	Negative	59	58	98
<b>TOTAL</b>		<b>120</b>	<b>119</b>	<b>99</b>

**Cocain***800 µg Benzoyllecgonin*

	Target value	Participant	Conform	Conform%
Architect	Positive	1	1	100
Cobas	Positive	4	4	100
Triage 8	Positive	8	8	100
Syva Rapid Test	Positive	7	7	100
Intex	Positive	8	8	100
Syntron Quik	Positive	1	1	100
Triage TOX	Positive	12	12	100
Other quantitativ	Positive	15	15	100
Other Rapid tests	Positive	60	56	93
<b>TOTAL</b>		<b>116</b>	<b>112</b>	<b>97</b>

**Methadon***1000 µg Methadon, 1000 µg EDDP*

	Target value	Participant	Conform	Conform%
Architect	Positive	1	1	100
Cobas	Positive	4	4	100
Triage 8	Positive	8	8	100
Intex	Positive	10	10	100
Syntron Quik	Positive	1	1	100
Triage TOX	Positive	12	12	100
Other quantitativ	Positive	14	14	100
Other Rapid tests	Positive	47	45	96
<b>TOTAL</b>		<b>97</b>	<b>95</b>	<b>98</b>

**EDDP***1000 µg EDDP*

	Target value	Participant	Conform	Conform%
All	Positive	1	1	100

**Methaqualon**

	Target value	Participant	Conform	Conform%
All	Negative	1	1	100

**Opiate***800 µg 6-Acetylmorphin*

	Target value	Participant	Conform	Conform%
Architect	Positive	1	1	100
Cobas	Positive	4	4	100
Triage 8	Positive	8	8	100
Syva Rapid Test	Positive	7	7	100
Intex	Positive	11	11	100
Syntron Quik	Positive	1	1	100
Triage TOX	Positive	12	12	100
Other quantitativ	Positive	16	16	100
Other Rapid tests	Positive	58	55	95
<b>TOTAL</b>		<b>118</b>	<b>115</b>	<b>97</b>

**6-Acetylmorphin***800 µg 6-Acetylmorphin*

	Target value	Participant	Conform	Conform%
All	Positive	2	2	100

**Paracetamol**

	Target value	Participant	Conform	Conform%
Triage TOX	Negative	11	11	100

**Phencyclidin**

	Target value	Participant	Conform	Conform%
Triage TOX	Negative	12	12	100
Other Rapid tests	Negative	8	8	100
<b>TOTAL</b>		<b>20</b>	<b>20</b>	<b>100</b>

**Tricyclic Antidepressiva**

	Target value	Participant	Conform	Conform%
Triage 8	Negative	12	12	100
Triage TOX	Negative	14	14	100
All methods	Negative	33	33	100
<b>TOTAL</b>		<b>12</b>	<b>12</b>	<b>100</b>

**LSD**

	Target value	Participant	Conform	Conform%
All	Negative	1	1	100

Zürich, 3.12.2016  
Dr. R. Fried