

U2 Urin test strip

MQ 2014/4

How to read the evaluation:

- The numbers in the boxes show how many participants have found a result.
- The box with a star is the target-value
- Your own result is in the box with a frame.

Method Combur Roche, visuell

Leucocytes	304*	120	6						
Lc/ul	neg	25	75		500				
Nitrite	5	434*							
	neg	pos							
pH	35	16	374*	11	2				
	5	5.5	6	6.5	7	7.5	8	8.5	9
Protein	7		339*	92	1				
g/L	neg		0.3	1	5				
Glucose	5	5	4	37	2	386*			
mmol/L	norm	2.8	5.5	17	28	55			
Ketone	426*		2	1		1			
	neg		+	++		+++			
Urobilinogene	405*	1							
μmol/L	norm	17		70	140	200			
Bilirubine	401*	3							
	neg	+	++	+++					
Hemoglobin	4	16	159*	154	99				
	neg	10	25	50	250				
Spec. gravity	2	1	21	202*	45	2			
	1000	1005	1010	1015	1020	1025	1030		

U2 Urin test strip

MQ 2014/4

How to read the evaluation:

- The numbers in the boxes show how many participants have found a result.
- The box with a star is the target-value
- Your own result is in the box with a frame.

Method Combur Urilux/Uriscys

Leucocytes	316*	23	1						
Lc/ul	neg	25	75		500				
Nitrite	1	339*							
	neg	pos							
pH	5	1	268*	62	1				
	5	5.5	6	6.5	7	7.5	8	8.5	9
Protein	10		221*	102	3	2			
g/L	neg		0.3	1	5				
Glucose	1	11	8	9	4	308*			
mmol/L	norm	2.8	5.5	17	28	55			
Ketone	336*		1						
	neg		+	++		+++			
Urobilinogene	337*	1							
μmol/L	norm	17		70	140	200			
Bilirubine	337*	1							
	neg	+	++	+++					
Hemoglobin	2	10	11	293*	23	1			
	neg	10	25	50	250				
Spec. gravity	2		85	203*	38	2			
	1000	1005	1010	1015	1020	1025	1030		

U2 Urin test strip

MQ 2014/4

How to read the evaluation

- The numbers in the boxes show how many participants have found a result.
- The box with a star is the target-value
- Your own result is in the box with a frame.

Method

Multistix Clinitek

Leucocytes	80*	5	1	1					
Lc/ul	neg	15	70	125	500				
Nitrite	1	86*							
	neg	pos							
pH			38	46*					
	5	5.5	6	6.5	7	7.5	8	8.5	9
Protein	1	1	13	71*	1				
g/L	neg	traccia	0.3	1	3	20			
Glucose	1			8	56*	20	2		
mmol/L	norm		5.5	14	28	55	111		
Ketone	85*								
	neg		+	++	+++	++++			
Urobilinogene	82*								
μmol/L	3.2	16	33	66	131				
Bilirubine	80*	1							
	neg	+	++	+++					
Hemoglobin		20	59*	7	1				
	neg	10	25	80	200				
Spec. gravity	1		3	21	52*	8			
	1000	1005	1010	1015	1020	1025	1030		

U2 Urin test strip

MQ 2014/4

How to read the evaluation:

- The numbers in the boxes show how many participants have found a result.
- The box with a star is the target-value
- Your own result is in the box with a frame.

Method Aution, PocketChem

Leucocytes	65*	1	1	1					
Lc/ul	neg	25	75		500				
Nitrite	2	66*							
	neg	pos							
pH			50*	18					
	5	5.5	6	6.5	7	7.5	8	8.5	9
Protein			39*	27	1	1			
g/L	neg		0.3	1	5				
Glucose			3	14	46*	5			
mmol/L	norm	2.8	5.5	17	28	55			
Ketone	66*			2					
	neg		+	++		+++			
Urobilinogene	68*								
μmol/L	norm	17		70	140	200			
Bilirubine	67*								
	neg	+	++	+++					
Hemoglobin	22	43*	2		1				
	neg	10	25	50	250				
Spec. gravity				40*	18	1			
	1000	1005	1010	1015	1020	1025	1030		

U2 Urin test strip

MQ 2014/4

How to read the evaluation:

- The numbers in the boxes show how many participants have found a result.
- The box with a star is the target-value
- Your own result is in the box with a frame.

Method

Combi Screen

Leucocytes	45*	1							
Lc/ul	neg	25	75		500				
Nitrite		46*							
	neg	pos							
pH	2		38*	2					
	5	5.5	6	6.5	7	7.5	8	8.5	9
Protein	2	4	35*	4	1				
g/L	neg		0.3	1	5				
Glucose	1	2		12	22*	9			
mmol/L	norm	2.8	5.5	17	28	55			
Ketone	41*					1			
	neg		+	++		+++			
Urobilinogene	40*								
μmol/L	norm	17		70	140	200			
Bilirubine	40*	1							
	neg	+	++	+++					
Hemoglobin	1	36*		5	3				
	neg	10	25	50	250				
Spec. gravity	1	1	4	23*	8	2			
	1000	1005	1010	1015	1020	1025	1030		