



Survey report MQ 2015-1 U3 Drugs of abuse

Amphetamine

	Target value	Participant	Conform	Conform%
Roche Cobas	Negative	9	9	100
Triage 8	Negative	9	9	100
Syva Rapid Test	Negative	4	4	100
Intex	Negative	6	6	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	11	11	100
Other quantitativ	Negative	12	11	92
Other Rapid tests	Negative	49	48	98
TOTAL		101	99	98

Methamphetamin

	Target value	Participant	Conform	Conform%
Intex	Negative	2	2	100
Syva Rapid Test	Negative	2	2	100
Triage TOX	Negative	11	11	100
Other Rapid tests	Negative	21	21	100
TOTAL		36	36	100

Ecstasy pos/neg

	Target value	Participant	Conform	Conform%
Other Rapid tests	Negative	8	8	100

Barbiturate

1000 µg Secobarbital / liter

	Target value	Participant	Conform	Conform%
Cobas	Positive	8	8	100
Triage 8	Positive	9	9	100
Intex	Positive	1	1	100
Triage TOX	Positive	11	11	100
Other quantitativ	Positive	8	8	100
Other Rapid tests	Positive	21	19	90
TOTAL		58	56	97

Benzodiazepine

	Target value	Participant	Conform	Conform%
Cobas	Negative	10	10	100
Triage 8	Negative	9	9	100
Intex	Negative	12	12	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	11	11	100
Other quantitativ	Negative	14	13	93
Other Rapid tests	Negative	56	56	100
TOTAL		113	112	99

Cannabinoide

	Target value	Participant	Conform	Conform%
Cobas	Negative	10	10	100
Triage 8	Negative	9	9	100
Syva Rapid Test	Negative	6	6	100
Intex	Negative	12	12	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	11	11	100
Other quantitativ	Negative	14	13	93
Other Rapid tests	Negative	60	60	100
TOTAL		123	122	99

Cocain

1000 µg Bencoylecgonin / liter

	Target value	Participant	Conform	Conform%
Cobas	Positive	10	10	100
Triage 8	Positive	9	9	100
Syva Rapid Test	Positive	6	6	100
Intex	Positive	12	12	100
Syntron Quik	Positive	1	1	100
Triage TOX	Positive	11	11	100
Other quantitativ	Positive	15	13	87
Other Rapid tests	Positive	61	61	100
TOTAL		125	123	98

LSD

5 µg LSD / liter

	Target value	Participant	Conform	Conform%
Cobas	Positive	1	1	100

Methadon

	Target value	Participant	Conform	Conform%
Cobas	Negative	10	10	100
Triage 8	Negative	9	9	100
Intex	Negative	13	13	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	11	11	100
Other quantitativ	Negative	14	12	86
Other Rapid tests	Negative	45	45	100
TOTAL		103	101	98

Opiate

	Target value	Participant	Conform	Conform%
Cobas	Negative	10	10	100
Triage 8	Negative	9	9	100
Syva Rapid Test	Negative	6	6	100
Intex	Negative	14	14	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	11	11	100
Other quantitativ	Negative	15	14	93
Other Rapid tests	Negative	61	61	100
TOTAL		127	126	99

Paracetamol

	Target value	Participant	Conform	Conform%
Triage TOX	Negative	11	11	100

Phencyclidin

	Target value	Participant	Conform	Conform%
Triage TOX	Negative	11	11	100
Other Rapid tests	Negative	8	8	100
TOTAL		15	15	100

Tricyclic Antidepressiva

	Target value	Participant	Conform	Conform%
Triage 8	Negative	8	8	100
Triage TOX	Negative	11	11	100
All methods	Negative	14	14	100
TOTAL		33	33	100

Zürich, 30.3.2015
Dr. R. Fried