



## Survey report MQ 2015-2 U3 Drugs of abuse

### Amphetamine

2000 µg D,L-Amphetamin / litre

	Target value	Participant	Conform	Conform%
Roche Cobas	Positive	8	7	88
Triage 8	Positive	9	8	89
Syva Rapid Test	Positive	4	3	75
Intex	Positive	5	5	100
Syntron Quik	Positive	1	1	100
Triage TOX	Positive	11	11	100
Other quantitativ	Positive	12	11	92
Other Rapid tests	Positive	51	44	86
<b>TOTAL</b>		<b>101</b>	<b>90</b>	<b>89</b>

### Methamphetamin

2000 µg D,L-Methamphetamin pro Liter

	Target value	Participant	Conform	Conform%
Intex	Positive	2	2	100
Syva Rapid Test	Positive	2	2	100
Triage TOX	Positive	11	4	36
Other Rapid tests	Positive	23	20	87
<b>TOTAL</b>		<b>38</b>	<b>28</b>	<b>74</b>

### Ecstasy pos/neg

	Target value	Participant	Conform	Conform%
Other Rapid tests	Negative	<b>10</b>	<b>10</b>	<b>100</b>

### Barbiturate

	Target value	Participant	Conform	Conform%
Cobas	Negativ	8	7	88
Triage 8	Negativ	9	8	89
Intex	Negativ	1	1	100
Triage TOX	Negativ	11	11	100
Other quantitativ	Negativ	8	8	100
Other Rapid tests	Negativ	24	23	96
<b>TOTAL</b>		<b>61</b>	<b>58</b>	<b>95</b>

**Benzodiazepine**

	Target value	Participant	Conform	Conform%
Cobas	Negative	9	8	89
Triage 8	Negative	9	8	89
Intex	Negative	11	11	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	11	11	100
Other quantitativ	Negative	14	14	100
Other Rapid tests	Negative	57	56	98
<b>TOTAL</b>		<b>112</b>	<b>109</b>	<b>97</b>

**Cannabinoide**

	Target value	Participant	Conform	Conform%
Cobas	Negative	9	8	89
Triage 8	Negative	9	8	89
Syva Rapid Test	Negative	6	6	100
Intex	Negative	11	11	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	11	11	100
Other quantitativ	Negative	15	15	100
Other Rapid tests	Negative	63	62	98
<b>TOTAL</b>		<b>125</b>	<b>122</b>	<b>98</b>

**Cocain**

	Target value	Participant	Conform	Conform%
Cobas	Negative	9	8	89
Triage 8	Negative	9	8	89
Syva Rapid Test	Negative	6	6	100
Intex	Negative	10	10	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	11	11	100
Other quantitativ	Negative	15	15	100
Other Rapid tests	Negative	64	63	98
<b>TOTAL</b>		<b>125</b>	<b>122</b>	<b>98</b>

**LSD**

	Target value	Participant	Conform	Conform%
Cobas	Negative	1	1	100

**Methadon**

	Target value	Participant	Conform	Conform%
Cobas	Negative	9	8	89
Triage 8	Negative	9	8	89
Intex	Negative	11	11	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	11	11	100
Other quantitativ	Negative	14	14	100
Other Rapid tests	Negative	48	44	92
<b>TOTAL</b>		<b>103</b>	<b>97</b>	<b>94</b>

**Opiate***1000 µg Morphin and 1000 µg Morphin-3-glucoronid (M3G) / litre*

	Target value	Participant	Conform	Conform%
Cobas	Positive	9	8	89
Triage 8	Positive	9	8	89
Syva Rapid Test	Positive	6	6	100
Intex	Positive	12	11	92
Syntron Quik	Positive	1	1	100
Triage TOX	Positive	11	11	100
Other quantitativ	Positive	15	15	100
Other Rapid tests	Positive	63	59	94
<b>TOTAL</b>		<b>126</b>	<b>119</b>	<b>94</b>

**Paracetamol**

	Target value	Participant	Conform	Conform%
Triage TOX	Negative	11	10	91

**Phencyclidin**

	Target value	Participant	Conform	Conform%
Triage TOX	Negative	11	11	100
Other Rapid tests	Negative	9	9	100
<b>TOTAL</b>		<b>20</b>	<b>20</b>	<b>100</b>

**Tricyclic Antidepressiva**

	Target value	Participant	Conform	Conform%
Triage 8	Negative	8	8	100
Triage TOX	Negative	11	11	100
All methods	Negative	15	15	100
<b>TOTAL</b>		<b>34</b>	<b>34</b>	<b>100</b>

Zürich, 6.7.2015  
 Dr. R. Fried