



## Survey report MQ 2015-3 U3 Drugs of abuse

### Amphetamine

	Target value	Participant	Conform	Conform%
Roche Cobas	Negativ	6	6	100
Triage 8	Negativ	9	9	100
Syva Rapid Test	Negativ	4	4	100
Intex	Negativ	5	5	100
Syntron Quik	Negativ	1	1	100
Triage TOX	Negativ	11	11	100
Other quantitativ	Negativ	11	11	100
Other Rapid tests	Negativ	50	50	100
<b>TOTAL</b>		<b>97</b>	<b>97</b>	<b>100</b>

### Methamphetamin

	Target value	Participant	Conform	Conform%
Intex	Negativ	2	2	100
Syva Rapid Test	Negativ	2	2	100
Triage TOX	Negativ	11	11	100
Other Rapid tests	Negativ	23	23	100
<b>TOTAL</b>		<b>38</b>	<b>38</b>	<b>100</b>

### Ecstasy pos/neg

	Target value	Participant	Conform	Conform%
Other Rapid tests	Negative	12	12	100

### Barbiturate

	Target value	Participant	Conform	Conform%
Cobas	Negativ	6	6	100
Triage 8	Negativ	9	9	100
Intex	Negativ	1	1	100
Triage TOX	Negativ	11	11	100
Other quantitativ	Negativ	8	8	100
Other Rapid tests	Negativ	24	24	100
<b>TOTAL</b>		<b>59</b>	<b>59</b>	<b>100</b>

**Benzodiazepine**

	Target value	Participant	Conform	Conform%
Cobas	Negative	7	7	100
Triage 8	Negative	9	9	100
Intex	Negative	11	11	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	11	11	100
Other quantitativ	Negative	13	12	92
Other Rapid tests	Negative	58	58	100
<b>TOTAL</b>		<b>110</b>	<b>109</b>	<b>99</b>

**Cannabinoide**

	Target value	Participant	Conform	Conform%
Cobas	Negative	8	8	100
Triage 8	Negative	9	9	100
Syva Rapid Test	Negative	6	6	100
Intex	Negative	12	11	92
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	11	11	100
Other quantitativ	Negative	14	14	100
Other Rapid tests	Negative	62	61	98
<b>TOTAL</b>		<b>123</b>	<b>121</b>	<b>98</b>

**Cocain**

	Target value	Participant	Conform	Conform%
Cobas	Negative	7	7	100
Triage 8	Negative	9	9	100
Syva Rapid Test	Negative	6	6	100
Intex	Negative	10	10	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	11	11	100
Other quantitativ	Negative	14	14	100
Other Rapid tests	Negative	63	63	100
<b>TOTAL</b>		<b>121</b>	<b>121</b>	<b>100</b>

**LSD**

	Target value	Participant	Conform	Conform%
Cobas	Negative	1	1	100

**Methadon***1500 µg/L Methadon, 750 µg/L EDDP*

	Target value	Participant	Conform	Conform%
Cobas	Positive	7	7	100
Triage 8	Positive	9	9	100
Intex	Positive	11	11	100
Syntron Quik	Positive	1	1	100
Triage TOX	Positive	11	11	100
Other quantitativ	Positive	13	13	100
Other Rapid tests	Positive	48	47	98
<b>TOTAL</b>		<b>100</b>	<b>99</b>	<b>99</b>

**Opiate**

	Target value	Participant	Conform	Conform%
Cobas	Negative	7	7	100
Triage 8	Negative	9	9	100
Syva Rapid Test	Negative	6	6	100
Intex	Negative	12	12	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	11	11	100
Other quantitativ	Negative	15	15	100
Other Rapid tests	Negative	61	61	100
<b>TOTAL</b>		<b>122</b>	<b>122</b>	<b>100</b>

**Paracetamol**

	Target value	Participant	Conform	Conform%
Triage TOX	Negative	11	11	100

**Phencyclidin**

	Target value	Participant	Conform	Conform%
Triage TOX	Negative	11	11	100
Other Rapid tests	Negative	9	9	100
<b>TOTAL</b>		<b>20</b>	<b>20</b>	<b>100</b>

**Tricyclic Antidepressiva***3000 µg/L Amitriptylin*

	Target value	Participant	Conform	Conform%
Triage 8	Positive	8	7	88
Triage TOX	Positive	11	11	100
All methods	Positive	16	16	100
<b>TOTAL</b>		<b>35</b>	<b>34</b>	<b>97</b>

Zürich, 28.9.2015

Dr. R. Fried