



Survey report MQ 2016-3 U3 Drugs of abuse

Amphetamine

	Target value	Participant	Conform	Conform%
Architect	Negative	1	1	100
Roche Cobas	Negative	3	3	100
Triage 8	Negative	8	8	100
Syva Rapid Test	Negative	2	2	100
Intex	Negative	6	6	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	12	12	100
Other quantitativ	Negative	12	12	100
Other Rapid tests	Negative	50	50	100
TOTAL		95	95	100

Methamphetamine

	Target value	Participant	Conform	Conform%
Intex	Negative	2	2	100
Syva Rapid Test	Negative	3	3	100
Triage TOX	Negative	12	12	100
Other Rapid tests	Negative	24	24	100
TOTAL		41	41	100

Ecstasy pos/neg

	Target value	Participant	Conform	Conform%
Other Rapid tests	Negative	12	12	100

Barbiturate

	Target value	Participant	Conform	Conform%
Architect	Negative	1	1	100
Cobas	Negative	3	3	100
Triage 8	Negative	8	8	100
Intex	Negative	1	1	100
Triage TOX	Negative	12	12	100
Other quantitativ	Negative	8	8	100
Other Rapid tests	Negative	24	24	100
TOTAL		57	57	100

Benzodiazepine

	Target value	Participant	Conform	Conform%
Architect	Negative	1	1	100
Cobas	Negative	4	4	100
Triage 8	Negative	8	8	100
Intex	Negative	11	11	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	12	12	100
Other quantitativ	Negative	13	13	100
Other Rapid tests	Negative	57	57	100
TOTAL		107	107	100

Cannabinoide

	Target value	Participant	Conform	Conform%
Architect	Negative	1	1	100
Cobas	Negative	6	6	100
Triage 8	Negative	8	8	100
Syva Rapid Test	Negative	6	6	100
Intex	Negative	12	12	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	12	12	100
Other quantitativ	Negative	14	14	100
Other Rapid tests	Negative	60	60	100
TOTAL		120	120	100

Cocain

	Target value	Participant	Conform	Conform%
Architect	Negative	1	1	100
Cobas	Negative	4	4	100
Triage 8	Negative	8	8	100
Syva Rapid Test	Negative	6	6	100
Intex	Negative	9	9	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	12	12	100
Other quantitativ	Negative	14	14	100
Other Rapid tests	Negative	61	60	98
TOTAL		116	115	99

Methadon

	Target value	Participant	Conform	Conform%
Architect	Negative	1	1	100
Cobas	Negative	4	4	100
Triage 8	Negative	8	8	100
Intex	Negative	11	11	100
Syntron Quik	Negative	1	1	100
Triage TOX	Negative	12	12	100
Other quantitativ	Negative	14	14	100
Other Rapid tests	Negative	49	49	100
TOTAL		100	100	100

Methaqualon*1000 µg Methaqualone*

	Target value	Participant	Conform	Conform%
All	Positive	1	1	100

Opiate*1000 µg Morphine*

	Target value	Participant	Conform	Conform%
Architect	Positive	1	1	100
Cobas	Positive	4	4	100
Triage 8	Positive	8	8	100
Syva Rapid Test	Positive	6	6	100
Intex	Positive	12	11	92
Syntron Quik	Positive	1	1	100
Triage TOX	Positive	12	12	100
Other quantitativ	Positive	15	15	100
Other Rapid tests	Positive	59	58	98
TOTAL		118	116	98

Paracetamol

	Target value	Participant	Conform	Conform%
Triage TOX	Negative	11	11	100

Phencyclidin*200 µg Phencyclidine*

	Target value	Participant	Conform	Conform%
Triage TOX	Positive	12	12	100
Other Rapid tests	Positive	8	8	100
TOTAL		20	20	100

Tricyclic Antidepressiva

	Target value	Participant	Conform	Conform%
Triage 8	Negative	7	7	100
Triage TOX	Negative	12	12	100
All methods	Negative	15	15	100
TOTAL		34	34	100

Zürich, 28.9.2016
Dr. R. Fried