

Report U03 Drugs of Abuse

Amphetamine

	Target Value	Participants	Conform	% Conform
Architect	positive	1	1	100.00
Cobas	positive	3	3	100.00
Syva Rapid Test	positive	2	1	50.00
Intex	positive	5	4	80.00
Syntron Quik	positive	1	1	100.00
Other quantitativ	positive	14	13	93.00
Others	positive	52	40	77.00
Total		78	63	81.00

Methamphetamine

	Target Value	Participants	Conform	% Conform
Intex	positive	2	2	100.00
Syva Rapid Test	positive	4	2	50.00
Others	positive	24	18	75.00
Total		30	22	73.00

Ecstasy pos/neg

	Target Value	Participants	Conform	% Conform
Others	negative	15	14	93.00

Barbiturates

	Target Value	Participants	Conform	% Conform
Architect	negative	1	1	100.00
Cobas	negative	3	3	100.00
Intex	negative	1	1	100.00
Other quantitativ	negative	10	10	100.00
Others	negative	25	25	100.00
Total		40	40	100.00

Benzodiazepines

	Target Value	Participants	Conform	% Conform
Architect	negative	1	1	100.00
Cobas	negative	4	4	100.00
Syva Rapid Test	negative	1	1	100.00
Intex	negative	10	10	100.00
Syntron Quik	negative	1	1	100.00
Other quantitativ	negative	15	15	100.00
Others	negative	58	57	98.00
Total		90	89	99.00

Cannabinoides

	Target Value	Participants	Conform	% Conform
Architect	positive	1	1	100.00
Cobas	positive	5	5	100.00
Syva Rapid Test	positive	7	6	86.00
Intex	positive	11	11	100.00
Syntron Quik	positive	1	1	100.00
Other quantitativ	positive	17	17	100.00
Others	positive	63	57	90.00
Total		105	98	93.00

Cocaine

	Target Value	Participants	Conform	% Conform
Architect	negative	1	1	100.00
Cobas	negative	4	4	100.00
Syva Rapid Test	negative	6	6	100.00
Intex	negative	8	8	100.00
Syntron Quik	negative	1	1	100.00
Other quantitativ	negative	16	16	100.00
Others	negative	62	61	98.00
Total		98	97	99.00

Methadone

	Target Value	Participants	Conform	% Conform
Architect	negative	1	1	100.00
Cobas	negative	4	4	100.00
Intex	negative	10	10	100.00
Syntron Quik	negative	1	1	100.00
Other quantitativ	negative	15	15	100.00
Others	negative	49	47	96.00
Total		80	78	98.00

Methaqualon

	Target Value	Participants	Conform	% Conform
all Participants	negative	1	1	100.00

EDDP

	Target Value	Participants	Conform	% Conform
all Participants	negative	1	1	100.00

Opiate

	Target Value	Participants	Conform	% Conform
Architect	negative	1	1	100.00
Cobas	negative	4	4	100.00
Syva Rapid Test	negative	6	6	100.00
Intex	negative	11	11	100.00
Syntron Quik	negative	1	1	100.00
Other quantitativ	negative	16	16	100.00
Others	negative	60	59	98.00
Total		99	98	99.00

Buprenorphin

	Target Value	Participants	Conform	% Conform
all Participants	negative	2	2	100.00

6-Acetylmorphin

	Target Value	Participants	Conform	% Conform
all Participants	negative	2	2	100.00

Phencyclidin

	Target Value	Participants	Conform	% Conform
Others	negative	7	7	100.00

TCA

	Target Value	Participants	Conform	% Conform
all Participants	negative	16	16	100.00

LSD

	Target Value	Participants	Conform	% Conform
Other quantitativ	negative	1	1	100.00

Commentary U3

To one liter of urine was added 350 ug 11-Nor-d9_carbonic acid, 350 ug d9-THC, 1500 ug D,L-methamphetamine and 1400 ug D,L-amphetamine. Correspondingly, cannabis, methamphetamine and amphetamine were supposed to be positive, all other analyses negative.