

Report U03 Drugs of Abuse

Amphetamine

	Target Value	Participants	Conform	% Conform
Architect	negative	1	1	100.00
Cobas	negative	3	3	100.00
Syva Rapid Test	negative	2	2	100.00
Intex	negative	5	5	100.00
Syntron Quik	negative	1	1	100.00
Other quantitativ	negative	14	14	100.00
Others	negative	52	51	98.00
Total		78	77	99.00

Methamphetamine

	Target Value	Participants	Conform	% Conform
Intex	negative	2	2	100.00
Syva Rapid Test	negative	4	4	100.00
Others	negative	22	21	95.00
Total		28	27	96.00

Ecstasy pos/neg

	Target Value	Participants	Conform	% Conform
Others	negative	13	12	92.00

Barbiturates

	Target Value	Participants	Conform	% Conform
Architect	negative	1	1	100.00
Cobas	negative	3	3	100.00
Intex	negative	1	1	100.00
Other quantitativ	negative	10	10	100.00
Others	negative	24	23	96.00
Total		39	38	97.00

Benzodiazepines

	Target Value	Participants	Conform	% Conform
Architect	positive	1	1	100.00
Cobas	positive	4	4	100.00
Syva Rapid Test	positive	1	1	100.00
Intex	positive	10	9	90.00
Syntron Quik	negative	1	1	100.00
Other quantitativ	positive	15	15	100.00
Others	positive	57	56	98.00
Total		89	87	98.00

Cannabinoïdes

	Target Value	Participants	Conform	% Conform
Architect	negative	1	1	100.00
Cobas	negative	4	4	100.00
Syva Rapid Test	negative	6	6	100.00
Intex	negative	11	11	100.00
Syntron Quik	negative	1	1	100.00
Other quantitativ	negative	17	17	100.00
Others	negative	62	61	98.00
Total		102	101	99.00

Cocaine

	Target Value	Participants	Conform	% Conform
Architect	negative	1	1	100.00
Cobas	negative	4	4	100.00
Syva Rapid Test	negative	5	5	100.00
Intex	negative	8	8	100.00
Syntron Quik	negative	1	1	100.00
Other quantitativ	negative	16	16	100.00
Others	negative	60	59	98.00
Total		95	94	99.00

Methadone

	Target Value	Participants	Conform	% Conform
Architect	positive	1	1	100.00
Cobas	positive	4	4	100.00
Intex	positive	10	9	90.00
Syntron Quik	positive	1	1	100.00
Other quantitativ	positive	15	15	100.00
Others	positive	48	46	96.00
Total		79	76	96.00

EDDP

	Target Value	Participants	Conform	% Conform
all Participants	positive	1	1	100.00

Opiate

	Target Value	Participants	Conform	% Conform
Architect	positive	1	1	100.00
Cobas	positive	4	4	100.00
Syva Rapid Test	positive	5	5	100.00
Intex	positive	11	9	82.00
Syntron Quik	positive	1	1	100.00
Other quantitativ	positive	16	16	100.00
Others	positive	59	57	97.00
Total		97	93	96.00

Buprenorphin

	Target Value	Participants	Conform	% Conform
all Participants	negative	2	2	100.00

6-Acetylmorphin

	Target Value	Participants	Conform	% Conform
all Participants	positive	2	2	100.00

Phencyclidin

	Target Value	Participants	Conform	% Conform
Others	negative	7	7	100.00

TCA

	Target Value	Participants	Conform	% Conform
all Participants	positive	15	14	93.00

LSD

	Target Value	Participants	Conform	% Conform
Other quantitativ	negative	1	1	100.00

Commentary U3

To a liter of urine are added 1000 µg Nordiazepam, 1000 µg Methadone, 650 µg EDDP, 2000 µg Nortriptyline and 900 µg 6-Acetylmorphine. Correspondingly, analyses for Benzodiazepines, Methadone, TCA and opiates should be positive, all others negative.