

U2 Bandelettes urinaires/HCG**Combur Roche, visuell**

Leucocytes	106*	400*	439*		9			
Lc/ μ L	neg	25	75		500			
Nitrite	922*	17						
	neg	pos						
pH	257*	145*	537*	7	3	1		
	5	5.5	6	6.5	7	7.5	8	8.5
Protéine	859*	1	84*	1				
g/L	neg		0.3	1	5			
Glucose	950*	1	1					
mmol/L	norm	2.8	5.5	17		55		
Cétone	939*	1	3			1		
mmol/L	neg		1	5		15		
Urobilinogène	881*							
μ mol/L	norm	17		70	140	200		
Bilirubine	877*	4						
	neg	+	++	+++	++++			
Erythrocytes/Hém	933*	9						
	neg	10	25	50	250			
Poids spécif.	7	4	31*	379*	167*	5	1	
	1000	1005	1010	1015	1020	1025	1030	

Commentaire U2

L'urine provenait de donneurs sains.

U2 Bandelettes urinaires/HCG**Combur Urilux/Urisys**

Leucocytes	55*	397*	503*	2	3			
Lc/ μ L	neg	25	100	250	500			
Nitrite	937*	20						
	neg	pos						
pH	153*	17*	685*	102*	1	1		
	5	5.5	6	6.5	7	7.5	8	8.5
Protéine	952*		7*	2				
g/L	neg		0.25	0.75	1.5	5		
Glucose	961*							
mmol/L	norm	3	6	17		56		
Cétone	940*	19	1	1				
mmol/L	neg	0.5	1.5	5		15		
Urobilinogène	955*	1						
μ mol/L	norm	17		70	140	200		
Bilirubine	958*							
	neg	+	++	+++				
Erythrocytes/Hém	946*	13	1	2				
Ery / μ L	neg	10	25	50	250			
Poids spécif.		3	145*	601*	182*	4	3	
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Commentaire U2

L'urine provenait de donneurs sains.

U2 Bandelettes urinaires/HCG**Multistix Clinitek**

Leucocytes	5*	110*	14*	4*				
Lc/μL	neg	15	70	125	500			
Nitrite	133*							
	neg	pos						
pH	3*	3*	113*	14*				
	5	5.5	6	6.5	7	7.5	8	8.5
Protéine	133*							
g/L	neg		Spur	0.3	1	3	20	
Glucose	133*							
mmol/L	neg		5.5	14	28	55	111	
Cétone	132*	1						
mmol/L	neg	0.5	1.5	4	8	16		
Urobilinogène	125*							
μmol/L	3.2	16	33	66	131			
Bilirubine	125*							
	neg	+	++	+++				
Erythrocytes/Hém	131*	2						
Ery / μl	neg	10	25	80	200			
Poids spécif.			7*	48*	66*	9*	1	
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Commentaire U2

L'urine provenait de donneurs sains.

U2 Bandelettes urinaires/HCG**Aution, PocketChem**

Leucocytes	10*	25*	121*	60*	6			
Lc/µL	neg	10-25	70-75	125-250	500			
Nitrite	222*							
	neg	pos						
pH	1*		204*	17*				
	5	5.5	6	6.5	7	7.5	8	8.5
Protéine	216*	6*						
g/L	neg	0.15	0.3	1	3-5	10-20		
Glucose	217*	4						
mmol/L	norm	2.2	5.6	14-17	28	55	111	
Cétone	219*	2		1				
mmol/L	neg	0.5	1.5	4-5	8-10	15-16		
Urobilinogène	221*							
µmol/L	norm	16	33-35	66-70	131-140	200		
Bilirubine	217*	3						
	neg	+	++	+++				
Erythrocytes/Hém	218*				1			
Ery / µl	neg	5-10	25	50-80	250			
Poids spécif.		4	1*	92*	95*	4*	3	
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Commentaire U2

L'urine provenait de donneurs sains.

U2 Bandelettes urinaires/HCG**Combi Screen**

Leucocytes	16*	34*	7*					
Lc/ μ L	neg	25	75			500		
Nitrite	58*							
	neg	pos						
pH	10*		46*	1*				
	5	5.5	6	6.5	7		8	
Protéine	55*	3*						
mg/dl	neg	Spur	30	100	500			
Glucose	57*							
mmol/L	neg	2.8	5.6	14	28	56		
Cétone	56*	1						
	neg	Spur	+	++	+++			
Urobilinogène	57*							
μ mol/L	norm	35	70	140	200			
Bilirubine	57*							
	neg	+	++	+++				
Erythrocytes/Hém	58*							
Ery / μ l	neg	5-10		50	300			
Poids spécif.			7*	42*	7*			
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Commentaire U2

L'urine provenait de donneurs sains.

U2 Bandelettes urinaires/HCG**Sysmex U**

Leucocytes	5*	66*	46*					
Lc/ μ L	neg	25	75		500			
Nitrite	116*							
	neg	pos						
pH	5*	32*	80*					
	5	5.5	6	6.5	7	7.5	8	8.5
Protéine	113*	3*						
g/l	neg	0.15	0.3	1.0	3.0	10.0		
Glucose	116*							
mmol/L	neg	2.8	5.6	14	28		111	
Cétone	116*							
mmol/L	neg		0.93	2.8	7.4			
Urobilinogène	115*							
μ mol/L	norm		34	68	135	202		
Bilirubine	115*							
	neg	+	++	+++				
Erythrocytes/Hém	112*							
Ery / μ l	neg	10	20	50	250			
Poids spécif.	2		53*	57*	4*			
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Commentaire U2

L'urine provenait de donneurs sains.