

U2 Urin Teststreifen/HCG**Combur Roche, visuell**

Leukozyten	882*	69*	1		1			
Lc/ μ L	neg	25	75		500			
Nitrit	937*	10						
	neg	pos						
pH	545*	202*	205*	3*	2		1	
	5	5.5	6	6.5	7	7.5	8	8.5
Protein	17		35*	320*	589*			
g/L	neg		0.3	1	5			
Glukose	945*	1	3	6		3		
mmol/L	norm	2.8	5.5	17		55		
Keton	943*	1	3			1		
mmol/L	neg		1	5		15		
Urobilinogen	892*	1						
μ mol/L	norm	17		70	140	200		
Bilirubin	877*	10	1					
	neg	+	++	+++	++++			
Erythrozyten/Häm	87*	678*	112*	59	21			
	neg	10	25	50	250			
spez. Gewicht	2	3	7*	62*	473*	52*	1	
	1000	1005	1010	1015	1020	1025	1030	

Kommentar U2

Der Urin enthielt Protein, Erythrozyten und HCG.

U2 Urin Teststreifen/HCG**Combur Urilux/Urissy**

Leukozyten	947*	13*	1		3			
Lc/ μ L	neg	25	100	250	500			
Nitrit	950*	17						
	neg	pos						
pH	509*	16*	427*	8*	1	1		
	5	5.5	6	6.5	7	7.5	8	8.5
Protein	1		11	38*	413*	505*		
g/L	neg		0.25	0.75	1.5	5		
Glukose	965*	2						
mmol/L	norm	3	6	17		56		
Keton	937*	24		3				
mmol/L	neg	0.5	1.5	5		15		
Urobilinogen	959*							
μ mol/L	norm	17		70	140	200		
Bilirubin	958*	2						
	neg	+	++	+++				
Erythrozyten/Häm	86*	567*	300*	12	2			
Ery / μ l	neg	10	25	50	250			
spez. Gewicht	1	2	13*	261*	609*	53*	4	
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin enthielt Protein, Erythrozyten und HCG.

U2 Urin Teststreifen/HCG

Multistix Clinitek

Leukozyten	133*	2*	1					
Lc/µL	neg	15	70	125	500			
Nitrit	135*	1						
	neg	pos						
pH	7*	9*	117*	2*			1	
	5	5.5	6	6.5	7	7.5	8	8.5
Protein	2		1	14*	94*	24*		
g/L	neg		Spur	0.3	1	3	20	
Glukose	136*							
mmol/L	neg		5.5	14	28	55	111	
Keton	135*	1						
mmol/L	neg	0.5	1.5	4	8	16		
Urobilinogen	125*	1						
µmol/L	3.2	16	33	66	131			
Bilirubin	129*							
	neg	+	++	+++				
Erythrozyten/Häm	8*	58*	68*	2				
Ery / µl	neg	10	25	80	200			
spez. Gewicht	1	1		8*	19*	82*	24*	
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin enthielt Protein, Erythrozyten und HCG.

U2 Urin Teststreifen/HCG**Aution, PocketChem**

Leukozyten	219*	6*		1	1			
Lc/ μ L	neg	10-25	70-75	125-250	500			
Nitrit	227*							
	neg	pos						
pH	2*	4*	179*	42*				
	5	5.5	6	6.5	7	7.5	8	8.5
Protein	2	2	9*	38*	165*	10*		
g/L	neg	0.15	0.3	1	3-5	10-20		
Glukose	223*	2			1			
mmol/L	norm	2.2	5.6	14-17	28	55	111	
Keton	221*	4		1				
mmol/L	neg	0.5	1.5	4-5	8-10	15-16		
Urobilinogen	226*		1					
μ mol/L	norm	16	33-35	66-70	131-140	200		
Bilirubin	222*	4						
	neg	+	++	+++				
Erythrozyten/Häm	56*	161*	8*					
Ery / μ l	neg	5-10	25	50-80	250			
spez. Gewicht	2		2	9*	145*	34*	14*	
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin enthielt Protein, Erythrozyten und HCG.

U2 Urin Teststreifen/HCG

Combi Screen

Leukozyten	56*	1*						
Lc/µL	neg	25	75			500		
Nitrit	57*							
	neg	pos						
pH	21*		34*					
	5	5.5	6	6.5	7		8	
Protein		1	1*	17*	37*			
mg/dl	neg	Spur	30	100	500			
Glukose	57*							
mmol/L	neg	2.8	5.6	14	28	56		
Keton	54*	1						
	neg	Spur	+	++	+++			
Urobilinogen	56*							
µmol/L	norm	35	70	140	200			
Bilirubin	56*							
	neg	+	++	+++				
Erythrozyten/Häm	54*	3*						
Ery / µl	neg	5-10		50	300			
spez. Gewicht	1			23*	27*	2*		
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin enthielt Protein, Erythrozyten und HCG.

U2 Urin Teststreifen/HCG

Sysmex U

Leukozyten	120*	6*	1					
Lc/µL	neg	25	75		500			
Nitrit	127*							
	neg	pos						
pH	6*	106*	16*					
	5	5.5	6	6.5	7	7.5	8	8.5
Protein	1		16*	8*	101*	1*		
g/l	neg	0.15	0.3	1.0	3.0	10.0		
Glukose	127*							
mmol/L	neg	2.8	5.6	14	28		111	
Keton	127*							
mmol/L	neg		0.93	2.8	7.4			
Urobilinogen	125*							
µmol/L	norm		34	68	135	202		
Bilirubin	125*	1						
	neg	+	++	+++				
Erythrozyten/Häm	5*	116*	2*					
Ery / µl	neg	10	20	50	250			
spez. Gewicht			73*	49*	5*			
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin enthielt Protein, Erythrozyten und HCG.