

**Report U03 Drugs of Abuse**

**Ethylglucuronid ql**

	Target Value	Participants	Conform	% Conform
Rapid Test	positive	2	2	100.00

**Amphetamine**

	Target Value	Participants	Conform	% Conform
Abbott	negative	4	4	100.00
Cobas	negative	5	5	100.00
Protzek	negative	4	4	100.00
Nal von Minden	negative	8	8	100.00
SureStep Abbott	negative	10	10	100.00
Intex	negative	3	3	100.00
Willi Fox	negative	4	4	100.00
MD Doctors Direct	negative	7	7	100.00
Clear Test	negative	5	5	100.00
Other quantitativ	negative	14	14	100.00
Others	negative	19	19	100.00
<b>Total</b>		<b>83</b>	<b>83</b>	<b>100.00</b>

**Methamphetamine**

	Target Value	Participants	Conform	% Conform
Syva Rapid Test	negative	2	2	100.00
Others	negative	24	24	100.00
<b>Total</b>		<b>26</b>	<b>26</b>	<b>100.00</b>

**Ecstasy pos/neg**

	Target Value	Participants	Conform	% Conform
Others	negative	20	20	100.00

**Barbiturates**

	Target Value	Participants	Conform	% Conform
Abbott	negative	3	3	100.00
Cobas	negative	5	5	100.00
Protzek	negative	2	2	100.00
Nal von Minden	negative	2	2	100.00
SureStep Abbott	negative	6	6	100.00
Willi Fox	negative	2	1	50.00
MD Doctors Direct	negative	3	3	100.00
Clear Test	negative	1	1	100.00
Other quantitativ	negative	7	7	100.00

Others	negative	12	12	100.00
<b>Total</b>		<b>43</b>	<b>42</b>	<b>98.00</b>

**Benzodiazepines**

	Target Value	Participants	Conform	% Conform
Abbott	positive	4	4	100.00
Cobas	positive	5	5	100.00
Protzek	positive	4	4	100.00
Nal von Minden	positive	8	7	88.00
SureStep Abbott	positive	10	6	60.00
Intex	positive	4	4	100.00
Willi Fox	positive	2	2	100.00
MD Doctors Direct	positive	8	3	38.00
Clear Test	positive	5	5	100.00
Other quantitativ	positive	16	16	100.00
Others	positive	17	11	65.00
<b>Total</b>		<b>83</b>	<b>67</b>	<b>81.00</b>

**Cannabinoïdes**

	Target Value	Participants	Conform	% Conform
Abbott	negative	5	5	100.00
Cobas	negative	5	5	100.00
Protzek	negative	4	4	100.00
Nal von Minden	negative	10	10	100.00
SureStep Abbott	negative	12	11	92.00
Syva Rapid Test	negative	2	2	100.00
Intex	negative	8	7	88.00
Willi Fox	negative	5	5	100.00
MD Doctors Direct	negative	9	9	100.00
Clear Test	negative	6	6	100.00
Other quantitativ	negative	15	15	100.00
Others	negative	26	26	100.00
<b>Total</b>		<b>107</b>	<b>105</b>	<b>98.00</b>

**Cocaine**

	Target Value	Participants	Conform	% Conform
Abbott	positive	4	4	100.00
Cobas	positive	5	5	100.00
Protzek	positive	4	4	100.00
Nal von Minden	positive	8	7	88.00
SureStep Abbott	positive	11	11	100.00
Syva Rapid Test	positive	1	1	100.00
Intex	positive	3	2	67.00
Willi Fox	positive	4	4	100.00
MD Doctors Direct	positive	8	8	100.00
Clear Test	positive	6	5	83.00
Other quantitativ	positive	16	15	94.00
Others	positive	21	21	100.00
<b>Total</b>		<b>91</b>	<b>87</b>	<b>96.00</b>

**Methadone**

	Target Value	Participants	Conform	% Conform
Abbott	negative	4	4	100.00
Cobas	negative	5	5	100.00
Protzek	negative	4	4	100.00
Nal von Minden	negative	6	6	100.00
SureStep Abbott	negative	9	9	100.00
Intex	negative	3	3	100.00
Willi Fox	negative	3	3	100.00
MD Doctors Direct	negative	9	9	100.00
Clear Test	negative	5	5	100.00
Other quantitativ	negative	12	12	100.00
Others	negative	20	19	95.00
<b>Total</b>		<b>80</b>	<b>79</b>	<b>99.00</b>

**Methaqualon**

	Target Value	Participants	Conform	% Conform
all Participants	negative	1	1	100.00

**EDDP**

	Target Value	Participants	Conform	% Conform
all Participants	negative	8	8	100.00

**Opiate**

	Target Value	Participants	Conform	% Conform
Abbott	positive	4	4	100.00
Cobas	positive	5	5	100.00
Protzek	positive	4	4	100.00
Nal von Minden	positive	8	7	88.00
SureStep Abbott	positive	9	9	100.00
Syva Rapid Test	positive	1	1	100.00
Intex	positive	3	3	100.00
Willi Fox	positive	4	4	100.00
MD Doctors Direct	positive	9	8	89.00
Clear Test	positive	6	6	100.00
Other quantitativ	positive	16	15	94.00
Others	positive	19	18	95.00
<b>Total</b>		<b>88</b>	<b>84</b>	<b>95.00</b>

**Buprenorphin**

	Target Value	Participants	Conform	% Conform
all Participants	negative	6	6	100.00

**6-Acetylmorphin**

	Target Value	Participants	Conform	% Conform
all Participants	positive	1	1	100.00

**Phencyclidin**

	Target Value	Participants	Conform	% Conform
Others	negative	4	4	100.00

**TCA**

	Target Value	Participants	Conform	% Conform
all Participants	positive	16	16	100.00

**LSD**

	Target Value	Participants	Conform	% Conform
Other quantitativ	negative	2	2	100.00

**Methamphetamine**

	Target Value	Participants	Conform	% Conform
Nal von Minden	negative	6	6	100.00

**Morphin**

	Target Value	Participants	Conform	% Conform
all Participants	positive	2	2	100.00

**Commentary U3**

*To one liter of urine 186 mg ethanol, 200 µg temazepam, 700 µg benzoylecgonine, 800 µg 6-acetylmorphine, 700 µg clomipramine, 900 µg desipramine, 400 µg amitriptyline, and 400 µg imipramine were added. Accordingly, benzodiazepines, opiates (incl. morphine, 6-acetylmorphine), cocaine and TCA should be positive, all other analyses negative. In Switzerland, the SCDAT recommendation is to report benzodiazepines at 100 µg as "positive". Individual tests unfortunately use higher limits. The sample also contained alcohol (4.8 mmol/l) and ethyl glucuronide (2.9 mg/l).*