

Report U03 Drugs of Abuse

Ethylglucuronid ql

	Target Value	Participants	Conform	% Conform
Rapid Test	positive	2	2	100.00

Amphetamine

	Target Value	Participants	Conform	% Conform
Abbott	negative	4	4	100.00
Cobas	negative	5	5	100.00
Protzek	negative	5	5	100.00
Nal von Minden	negative	8	8	100.00
SureStep Abbott	negative	12	12	100.00
Intex	negative	3	3	100.00
Willi Fox	negative	3	3	100.00
MD Doctors Direct	negative	7	7	100.00
Clear Test	negative	4	4	100.00
Other quantitativ	negative	12	12	100.00
Others	negative	20	18	90.00
Total		83	81	98.00

Methamphetamine

	Target Value	Participants	Conform	% Conform
Syva Rapid Test	negative	2	2	100.00
Others	negative	28	27	96.00
Total		30	29	97.00

Ecstasy pos/neg

	Target Value	Participants	Conform	% Conform
Others	negative	21	21	100.00

Barbiturates

	Target Value	Participants	Conform	% Conform
Abbott	negative	3	3	100.00
Cobas	negative	5	5	100.00
Protzek	negative	2	2	100.00
Nal von Minden	negative	2	2	100.00
SureStep Abbott	negative	8	8	100.00
Willi Fox	negative	1	1	100.00
MD Doctors Direct	negative	3	3	100.00
Clear Test	negative	1	1	100.00
Other quantitativ	negative	6	6	100.00

Others	negative	13	11	85.00
Total		44	42	95.00

Benzodiazepines

	Target Value	Participants	Conform	% Conform
Abbott	negative	4	4	100.00
Cobas	negative	5	5	100.00
Protzek	negative	4	4	100.00
Nal von Minden	negative	8	8	100.00
SureStep Abbott	negative	11	11	100.00
Intex	negative	4	4	100.00
Willi Fox	negative	2	2	100.00
MD Doctors Direct	negative	8	8	100.00
Clear Test	negative	4	4	100.00
Other quantitativ	negative	15	14	93.00
Others	negative	18	17	94.00
Total		83	81	98.00

Cannabinoide

	Target Value	Participants	Conform	% Conform
Abbott	positive	5	5	100.00
Cobas	positive	5	5	100.00
Protzek	positive	5	4	80.00
Nal von Minden	positive	11	9	82.00
SureStep Abbott	positive	14	11	79.00
Syva Rapid Test	positive	2	1	50.00
Intex	positive	7	5	71.00
Willi Fox	positive	4	4	100.00
MD Doctors Direct	positive	9	9	100.00
Clear Test	positive	5	5	100.00
Other quantitativ	positive	14	14	100.00
Others	positive	27	24	89.00
Total		108	96	89.00

Cocaine

	Target Value	Participants	Conform	% Conform
Abbott	negative	4	4	100.00
Cobas	negative	5	5	100.00
Protzek	negative	5	5	100.00
Nal von Minden	negative	8	6	75.00
SureStep Abbott	negative	13	13	100.00
Syva Rapid Test	negative	1	1	100.00
Intex	negative	3	3	100.00
Willi Fox	negative	3	3	100.00
MD Doctors Direct	negative	8	8	100.00
Clear Test	negative	5	5	100.00
Other quantitativ	negative	14	14	100.00
Others	negative	22	20	91.00
Total		91	87	96.00

Methadone

	Target Value	Participants	Conform	% Conform
Abbott	negative	4	4	100.00
Cobas	negative	5	5	100.00
Protzek	negative	4	4	100.00
Nal von Minden	negative	6	6	100.00
SureStep Abbott	negative	11	11	100.00
Intex	negative	3	3	100.00
Willi Fox	negative	2	2	100.00
MD Doctors Direct	negative	8	8	100.00
Clear Test	negative	4	4	100.00
Other quantitativ	negative	11	11	100.00
Others	negative	21	19	90.00
Total		79	77	97.00

Methaqualon

	Target Value	Participants	Conform	% Conform
all Participants	positive	1	1	100.00

EDDP

	Target Value	Participants	Conform	% Conform
all Participants	negative	8	7	88.00

Opiate

	Target Value	Participants	Conform	% Conform
Abbott	positive	4	4	100.00
Cobas	positive	5	5	100.00
Protzek	positive	5	4	80.00
Nal von Minden	positive	8	6	75.00
SureStep Abbott	positive	11	9	82.00
Syva Rapid Test	positive	1	1	100.00
Intex	positive	3	3	100.00
Willi Fox	positive	3	3	100.00
MD Doctors Direct	positive	8	8	100.00
Clear Test	positive	5	5	100.00
Other quantitativ	positive	14	14	100.00
Others	positive	20	18	90.00
Total		87	80	92.00

Buprenorphin

	Target Value	Participants	Conform	% Conform
all Participants	negative	7	6	86.00

6-Acetylmorphin

	Target Value	Participants	Conform	% Conform
all Participants	negative	1	1	100.00

Phencyclidin

	Target Value	Participants	Conform	% Conform
Others	negative	4	4	100.00

TCA

	Target Value	Participants	Conform	% Conform
all Participants	negative	18	17	94.00

LSD

	Target Value	Participants	Conform	% Conform
Other quantitativ	negative	2	2	100.00

Methamphetamine

	Target Value	Participants	Conform	% Conform
Nal von Minden	negative	6	6	100.00

Morphin

	Target Value	Participants	Conform	% Conform
all Participants	positive	3	3	100.00

Commentary U3

600 µg methaqualone, 800 µg morphine, 600 µg morphine-3-glucuronide, 200 µg THC-COOH, 200 µg CBD and approx. 1000 µg ethyl-?-D-glucuronide were added to one liter of urine.