

U2 Urin Teststreifen/HCG**Combur Roche, visuell**

Leukozyten	51	122*	578*		55*			
Lc/ μ L	neg	25	75		500			
Nitrit	7	803*						
	neg	pos						
pH	65	9*	671*	50*	8			
	5	5.5	6	6.5	7	7.5	8	8.5
Protein	565*		239*	5	1			
g/L	neg		0.3	1	5			
Glukose	813*	1						
mmol/L	norm	2.8	5.5	17		55		
Keton	803*		1					
mmol/L	neg		1	5		15		
Urobilinogen	746*							
μ mol/L	norm	17		70	140	200		
Bilirubin	742*	5						
	neg	+	++	+++	++++			
Erythrozyten/Häm	6	4	74*	334*	392*			
Ery / μ L	neg	10	25	50	250			
spez. Gewicht	4	12	181*	245*	12*	3	1	
	1000	1005	1010	1015	1020	1025	1030	

Kommentar U2

Der Urin stammte von gesunden Spendern, dem Blut und Nitrit hinzugefügt wurde.

U2 Urin Teststreifen/HCG**Combur Urilux/UriSYS**

Leukozyten	22	100*	946*		81*			
Lc/µL	neg	25	100		500			
Nitrit	3	1147*						
	neg	pos						
pH	8	4*	538*	521*	72*		1	
	5	5.5	6	6.5	7	7.5	8	8.5
Protein	1119*		25*	1		1		
g/L	neg		0.25	0.75	1.5	5		
Glukose	1142*	1	1					
mmol/L	norm	3	6	17		56		
Keton	1142*	2		1				
mmol/L	neg	0.5	1.5	5		15		
Urobilinogen	1141*							
µmol/L	norm	17		70	140	200		
Bilirubin	1140*	1						
	neg	+	++	+++				
Erythrozyten/Häm	3	2	34*	85*	1025*			
Ery / µl	neg	10	25	50	250			
spez. Gewicht		62*	727*	293*	15		1	
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin stammte von gesunden Spendern, dem Blut und Nitrit hinzugefügt wurde.

U2 Urin Teststreifen/HCG**Multistix Clinitek**

Leukozyten	2	22*	109*	6*	2			
Lc/ μ L	neg	15	70	125	500			
Nitrit		141*						
	neg	pos						
pH		1	26*	110*	3*			
	5	5.5	6	6.5	7	7.5	8	8.5
Protein	140*		1					
g/L	neg	Spur	0.3	1	3	20		
Glukose	140*							
mmol/L	neg		5.5	14	28	55	111	
Keton	140*							
mmol/L	neg	0.5	1.5	4	8	16		
Urobilinogen	132*	1						
μ mol/L	3.2	16	33	66	131			
Bilirubin	134*							
	neg	+	++	+++				
Erythrozyten/Häm	1	1	4*	120*	15*			
Ery / μ L	neg	10	25	80	200			
spez. Gewicht		1	32*	76*	24*	3		
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin stammte von gesunden Spendern, dem Blut und Nitrit hinzugefügt wurde.

U2 Urin Teststreifen/HCG**Aution, PocketChem**

Leukozyten	6	5	129*	108*	66*			
Lc/µL	neg	10-25	70-75	125-250	500			
Nitrit	2	311*						
	neg	pos						
pH			41*	270*				
	5	5.5	6	6.5	7	7.5	8	8.5
Protein	308*	2*	2					
g/L	neg	0.15	0.3	1	3-5	10-20		
Glukose	309*	1	1					
mmol/L	norm	2.8	5.6	14-17	28	55	111	
Keton	311*	1						
mmol/L	neg	0.5	1.5	4-5	8-10	15-16		
Urobilinogen	313*		1					
µmol/L	norm	16	33-35	66-70	131-140	200		
Bilirubin	306*	3	1					
	neg	+	++	+++				
Erythrozyten/Häm	11	6	6*	221*	66*			
Ery / µl	neg	5-10	25	50-80	250			
spez. Gewicht		1	48*	181*	37*	2		
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin stammte von gesunden Spendern, dem Blut und Nitrit hinzugefügt wurde.

U2 Urin Teststreifen/HCG**Sysmex U**

Leukozyten	2	14*	243*		2*			
Lc/μL	neg	25	75		500			
Nitrit	1	259*						
	neg	pos						
pH	4	16*	201*	39*	1*			
	5	5.5	6	6.5	7	7.5	8	8.5
Protein	259*		1					
g/l	neg	0.15	0.3	1.0	3.0	10.0		
Glukose	257*		1					
mmol/L	neg	2.8	5.6	14	28		111	
Keton	259*							
mmol/L	neg		0.93	2.8	7.4			
Urobilinogen	258*		1					
μmol/L	norm		34	68	135	202		
Bilirubin	258*	1						
	neg	+	++	+++				
Erythrozyten/Häm	7	1	5*	225*	20*			
Ery / μl	neg	10	20	50	250			
spez. Gewicht	1	9	200*	44*	2*			
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin stammte von gesunden Spendern, dem Blut und Nitrit hinzugefügt wurde.

U2 Urin Teststreifen/HCG**Combi Screen**

Leukozyten	3	14*	41*			3*		
Lc/µL	neg	25	75			500		
Nitrit	1	60*						
	neg	pos						
pH			58*	2*	1*			
	5	5.5	6	6.5	7		8	
Protein	59*	2*						
mg/dl	neg	Spur	30	100	500			
Glukose	61*							
mmol/L	neg	2.8	5.6	14	28	56		
Keton	59*	1	1					
	neg	Spur	+	++	+++			
Urobilinogen	61*							
µmol/L	norm	35	70	140	200			
Bilirubin	60*							
	neg	+	++	+++				
Erythrozyten/Häm		2*		53*	6*			
Ery / µl	neg	5-10		50	300			
spez. Gewicht		2	17*	36*	2*	1		
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin stammte von gesunden Spendern, dem Blut und Nitrit hinzugefügt wurde.

U2 Urin Teststreifen/HCG**Medi-Test/Urixxon**

Leukozyten		2*	10*					
Lc/µL	neg	10-25	70-75	125-250	500			
Nitrit		12*						
	neg	pos						
pH		1*	2*	8*	1*			
	5	5.5	6	6.5	7	7.5	8	8.5
Protein	11*		1*					
g/L	neg	0.15	0.3	1	3-5	10-20		
Glukose	12*							
mmol/L	neg	2.2	5.6	14-17	28	55	111	
Keton	12*							
mmol/L	neg	1.5	4-5	8-10	15-16			
Urobilinogen	13*							
µmol/L	norm	16	33-35	66-70	131-140	200		
Bilirubin	12*							
	neg	+	++	+++				
Erythrozyten/Häm			1*	6*	5*			
Ery / µl	neg	5-10	25	50-80	250			
spez. Gewicht		2*	8*					
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin stammte von gesunden Spendern, dem Blut und Nitrit hinzugefügt wurde.