

U2 Urin Teststreifen/HCG**Combur Roche, visuell**

Leukozyten	451*	303*	43					
Lc/μL	neg	25	75		500			
Nitrit	786*	11						
	neg	pos						
pH	81*	13*	650*	25*	23		2	
	5	5.5	6	6.5	7	7.5	8	8.5
Protein	8		154*	570*	67*			
g/L	neg		0.3	1	5			
Glukose	5	6	29*	509*		256*		
mmol/L	norm	2.8	5.5	17		55		
Keton	790*		3					
mmol/L	neg		1	5		15		
Urobilinogen	742*	2						
μmol/L	norm	17		70	140	200		
Bilirubin	741*	3						
	neg	+	++	+++	++++			
Erythrozyten/Häm	284*	443*	47*	18	7			
Ery / μl	neg	10	25	50	250			
spez. Gewicht	9	21*	141*	257*	18*	4		
	1000	1005	1010	1015	1020	1025	1030	

Kommentar U2

Der Urin stammte von gesunden Spendern, dem Glukose und Protein hinzugefügt wurde.

U2 Urin Teststreifen/HCG**Combur Urilux/UriSYS**

Leukozyten	933*	226*	12					
Lc/µL	neg	25	100		500			
Nitrit	1135*	38						
	neg	pos						
pH	8	2*	874*	239*	50*			
	5	5.5	6	6.5	7	7.5	8	8.5
Protein	3		24*	153*	977*	14*		
g/L	neg		0.25	0.75	1.5	5		
Glukose	4	2	28*	1017*		122*		
mmol/L	norm	3	6	17		56		
Keton	1167*	1	1					
mmol/L	neg	0.5	1.5	5		15		
Urobilinogen	1167*							
µmol/L	norm	17		70	140	200		
Bilirubin	1165*							
	neg	+	++	+++				
Erythrozyten/Häm	585*	483*	91*	13	1			
Ery / µl	neg	10	25	50	250			
spez. Gewicht	7	129*	670*	293*	20*		1	
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin stammte von gesunden Spendern, dem Glukose und Protein hinzugefügt wurde.

U2 Urin Teststreifen/HCG**Multistix Clinitek**

Leukozyten	104*	32*		2				
Lc/µL	neg	15	70	125	500			
Nitrit	138*							
	neg	pos						
pH		1*	48*	71*	14*			
	5	5.5	6	6.5	7	7.5	8	8.5
Protein			3*	87*	46*	1*		
g/L	neg	Spur	0.3	1	3	20		
Glukose			3	115*	15*	5*		
mmol/L	neg		5.5	14	28	55	111	
Keton	136*	1						
mmol/L	neg	0.5	1.5	4	8	16		
Urobilinogen	132*							
µmol/L	3.2	16	33	66	131			
Bilirubin	132*							
	neg	+	++	+++				
Erythrozyten/Häm	14*	114*	11*					
Ery / µl	neg	10	25	80	200			
spez. Gewicht			26*	72*	32*	3*		
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin stammte von gesunden Spendern, dem Glukose und Protein hinzugefügt wurde.

U2 Urin Teststreifen/HCG**Aution, PocketChem**

Leukozyten	271*	45*	6	2				
Lc/µL	neg	10-25	70-75	125-250	500			
Nitrit	322*							
	neg	pos						
pH			41*	281*	1*			
	5	5.5	6	6.5	7	7.5	8	8.5
Protein		1	32*	288*	3*			
g/L	neg	0.15	0.3	1	3-5	10-20		
Glukose		1	14*	141*	168*			
mmol/L	norm	2.8	5.6	14-17	28	55	111	
Keton	316*	6						
mmol/L	neg	0.5	1.5	4-5	8-10	15-16		
Urobilinogen	321*	3	1					
µmol/L	norm	16	33-35	66-70	131-140	200		
Bilirubin	316*	3						
	neg	+	++	+++				
Erythrozyten/Häm	175*	134*	9*	1				
Ery / µl	neg	5-10	25	50-80	250			
spez. Gewicht		4	57*	182*	32*	3*	1	
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin stammte von gesunden Spendern, dem Glukose und Protein hinzugefügt wurde.

U2 Urin Teststreifen/HCG**Sysmex U**

Leukozyten	206*	57*	9		1			
Lc/µL	neg	25	75		500			
Nitrit	273*							
	neg	pos						
pH	9*	27*	230*	8*				
	5	5.5	6	6.5	7	7.5	8	8.5
Protein		3	6*	265*	1*			
g/l	neg	0.15	0.3	1.0	3.0	10.0		
Glukose	4	1	7*	213*	48*		2	
mmol/L	neg	2.8	5.6	14	28		111	
Keton	271*							
mmol/L	neg		0.93	2.8	7.4			
Urobilinogen	272*		2					
µmol/L	norm		34	68	135	202		
Bilirubin	268*	3						
	neg	+	++	+++				
Erythrozyten/Häm	28*	223*	21*	1				
Ery / µl	neg	10	20	50	250			
spez. Gewicht	1*	50*	207*	9*	1		1	
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin stammte von gesunden Spendern, dem Glukose und Protein hinzugefügt wurde.

U2 Urin Teststreifen/HCG**Combi Screen**

Leukozyten	64*	4*						
Lc/µL	neg	25	75			500		
Nitrit	66*	2						
	neg	pos						
pH	1		64*	2*				
	5	5.5	6	6.5	7		8	
Protein		3	33*	30*	1*			
mg/dl	neg	Spur	30	100	500			
Glukose			3*	36*	28*	1		
mmol/L	neg	2.8	5.6	14	28	56		
Keton	66*	2						
	neg	Spur	+	++	+++			
Urobilinogen	67*							
µmol/L	norm	35	70	140	200			
Bilirubin	66*							
	neg	+	++	+++				
Erythrozyten/Häm	58*	9*		1*				
Ery / µl	neg	5-10		50	300			
spez. Gewicht		3*	41*	19*	1*			
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin stammte von gesunden Spendern, dem Glukose und Protein hinzugefügt wurde.

U2 Urin Teststreifen/HCG**Medi-Test/Urixxon**

Leukozyten	11*	1*						
Lc/µL	neg	10-25	70-75	125-250	500			
Nitrit	12*							
	neg	pos						
pH			10*	2*				
	5	5.5	6	6.5	7	7.5	8	8.5
Protein			3*	6*	3*			
g/L	neg	0.15	0.3	1	3-5	10-20		
Glukose			6*	3*	3*			
mmol/L	neg	2.2	5.6	14-17	28	55	111	
Keton	12*							
mmol/L	neg	1.5	4-5	8-10	15-16			
Urobilinogen	13*							
µmol/L	norm	16	33-35	66-70	131-140	200		
Bilirubin	12*							
	neg	+	++	+++				
Erythrozyten/Häm	4*	7*		1*				
Ery / µl	neg	5-10	25	50-80	250			
spez. Gewicht		2*	6*	2*				
	1.000	1.005	1.010	1.015	1.020	1.025	1.030	

Kommentar U2

Der Urin stammte von gesunden Spendern, dem Glukose und Protein hinzugefügt wurde.