

Methode

Aution Pocketchem

	1	2	3	4	5				
Leukozyten	Neg	10-25lc/ μ l 25	70-75lc/ μ l 75	125-250 250	500lc/ μ l 500				
	1	2							
Nitrit	Neg	Pos							
	1	2	3	4	5	6	7	8	9
pH	5	5.5	6	6.5	7	7.5	8	8.5	9
	1	2	3	4	5			6	
Protein	Neg	0.15g/L 13mg/dl Spur (+/-)	0.3g/L 30mg/dl +	1g/L 100mg/dl 2+	3-5g/L 300mg/dl 3+			10-20g/L 1000mg/dl 4+	
	1	2	3	4	5	6	7		
Glukose	norm	2.8 mmol/L 50 mg/dl Spur (+/-)	5.6 mmol/L 100 mg/dl +	11 mmol/L 200mg/dl 2+	28 mmol/L 500mg/dl 3+	55 mmol/L 1000mg/dl 4+	111 mmol/L 2000 5+		
	1	2	3			4	5	6	
Keton	Neg	0.5mM 5 mg/dl Spur (+/-)	1.5mM 15mg/dl 1+			4-5 mM 40mg/dl 2+	8-10 mM 80mg/dl 3+	15-16 mM 150mg/dl 4+	
	1	2	3			4	5	6	
Urobilinogen	Norm	16 μ mol/l 1 mg/dl +/-	33-35 μ M 2 mg/dl 1+			66-70 μ M 4 mg/dl 2+	140 μ M 8 mg/dl 3+	200 μ mol/l > 8 mg/dl 4+	
	1	2	3			4			
Bilirubin	Neg	+ 0.5mg/dl	2+ 2 mg/dl			3+ 6 mg/dl			
	1	2	3			4			
HB / Ery	Neg	5-10 Ec/ μ l 0.03 mg/dl +/-	25 Ec/ μ l 0.06 mg/dl +			50-80 Ec/ μ l 0.2 mg/dl 2+			
	1	2	3	4	5	6	7		
Spez. Gewicht	1.000	1.005	1.010	1.015	1.020	1.025	1.030		