

Methode

Siemens Multistix visuell oder clinitek

	1	2	3	4	5			
Leukozyten	Neg	15 Spur	70 +	125 2+	500 3+			
	1	2						
Nitrit	Neg	Pos						
	1	2	3	4	5	6	7	8
pH	5	5.5	6	6.5	7	7.5	8	8.5
	1		2	3	4	5	6	
Protein	Neg		Spur 15 mg/dl (+/-)	0.3g/L 30mg/dl +	1g/L 100mg/dl 2+	3g/L 300mg/dl 3+	≥20g/L ≥2000mg/dl 4+	
	1		3	4	5	6	7	
Glukose	Neg		5.5mmol/L 1g/L 100mg/dl Spur	14mmol/L 2.5g/L 250mg/dl +	28mmol/L 5g/L 500mg/dl 2+	55mmol/L 10g/L 1000mg/dl 3+	≥111mmol/L ≥20g/L ≥2000mg/dl 4+	
	1	2	3	4	5	6		
Keton	Neg	0.5mmol/L 0.05g/L 5mg/dl Spur	1.5mmol/L 0.15g/L 15mg/dl +	4mmol/L 0.4g/L 40mg/dl 2+	8mmol/L 0.8g/L 80mg/dl 3+	≥16mmol/L ≥1.6g/L ≥160mg/dl 4+		
	1	2	3	4	5			
Urobilinogen	3.2µmol/L 0.2mg/dl norm	16µmol/L 1mg/dl	33µmol/L 2mg/dl	66µmol/L 4mg/dl	≥131µmol/L ≥8mg/dl			
	1	2	3	4				
Bilirubin	Neg	1+	2+	3+				
	1	2	3	4	5			
HB / ERY	Neg	10Ery/ µL Spur	25Ery/µL +	80Ery/µL 2+	200Ery/µL 3+			
	1	2	3	4	5	6	7	
Spez. Gewicht	1.000	1.005	1.010	1.015	1.020	1.025	1.030	
	1	2						
HCG	Neg	Pos						

Methode: