

Methode

Roche, Combur visuell

Leukozyten	1 Neg	2 25lc/ μ l +	3 75lc/ μ l 2+				5 500lc/ μ l 3+		
Nitrit	1 Neg	2 Pos							
pH	1 5	2 5.5	3 6	4 6.5	5 7	6 7.5	7 8	8 8.5	9 9
Protein	1 Neg			3 0.3g/l 30mg/dl +	4 1g/l 100mg/dl 2+	5 5g/l 500mg/dl 3+			
Glukose	1 Normal		2 2.8mmol/L 50mg/dl +		3 5.5mmol/L 100mg/dl 2+		4 17mmol/L 300mg/dl 3+		6 55mmol/L 1000mg/dl 4+
Keton	1 Neg		3 1mmol/l 10mg/dl +		4 5mmol/L 50mg/dl 2+		6 15mmol/L 150mg/dl 3+		
Urobilinogen	1 Normal		2 17 μ mol/L 1mg/dl +		4 70 μ mol/L 4mg/dl 2+		5 140 μ mol/L 8mg/dl 3+		6 200 μ mol/L 12mg/dl 4+
Bilirubin	1 Neg		2 1+		3 2+		4 3+		
HB / Ery	1 Neg		2 5-10Ery/ μ L +		3 25Ery/ μ L 2+		4 50Ery/ μ L 3+		5 250Ery/ μ L 4+
Spez. Gewicht	1 1.000	2 1.005	3 1.010	4 1.015	5 1.020	6 1.025	7 1.030		